

Details



Simple. Fast. Accurate.



PRESCRIBING TOOLS

The m'eyeFit® mirror system has a full suite of sales support and patient education modules to support you and your team in the dispensary.

Product Demos

Lens Designs | Materials
No-Glare | Polarized
Photochromic | Anti-Fog

Frame Styling

Compare multiple frame styles quickly and easily

m'eyeFit® mirror specifications

Windows-based personal computer **required** (not included):

OS: Microsoft Windows XP, minimum;
Windows 7 or 10, recommended

Processor: Dual-core 1.8GHz, minimum
RAM 2GB DDR2, minimum

USB: One USB 2.0 port free or higher

Internet access for registration and updates recommended

Dimensions:

Mirror: 11.4" (H) x 4" (W)

Base of mirror: 7" (W) x 5.6" (D)

Mirror USB cord approximately 6'

m'eyeFit® mirror system



service & support

Your purchase includes:

- A limited three-year warranty
- A dedicated team of support staff
- Automatic software and service enhancements
- Access training videos by searching "meyeFitUSA" on YouTube.

Contact your Essilor Sales Consultant to purchase your m'eyeFit® mirror system today!

Technical support (877) 457-3760 or
meyefit@essilorusa.com



All it takes is a look in the mirror



Precision and Personalization

m'eyeFit® mirror system is a flexible solution designed to match different practice needs.

Using cutting edge 3D modeling technology to ensure accurate measurements, it is a compact plug and play solution that adapts to all dispensing tables. It's simple, fast and compatible with all Windows-based PCs.

CONFIGURATION

PC required for operation

- Simple plug and play into your existing PC*
- Includes software for installation on your PC

**Windows-based systems only.*



m'eyeFit software on PC (PC not included)

Why Choose a m'eyeFit® mirror digital measuring system?

• Simple. Fast. Accurate.

The m'eyeFit® mirror digital measuring system offers reliable solutions that deliver a modern patient experience. The plug and play designs make the device easy to set up for simple, fast, and accurate measurements.

• Reliable

With the m'eyeFit® mirror system, your patient records are always accessible. PC-driven processing capability internalizes the measurement calculations rather than sending them out to a remote server or into the "Cloud." This eliminates susceptibility to issues like server outages, lost internet connections and application failures.

BENEFITS

Patient

- Quick, comfortable process
- Precise measurements for the latest lens technology
- Personalized visual solutions for optimized vision

Eyecare Professional

- Simple and reliable measurement protocol
- Easy integration into current dispensing workflow
- Powerful processing capability
- Dedicated service and support
- Access to advanced personalized lenses for both single vision and progressive prescriptions

MEASUREMENTS

The m'eyeFit® mirror system accurately measures for a full range of lens designs, from single vision to personalized digital progressive designs.

Standard Fitting

Monocular PD
 Total PD
 Fitting Heights
 Segment Heights
 Bifocal Heights
 A, B, DBL

Frame Fitting Personalization

Vertex Distance
 Wrap Angle
 Pantoscopic Tilt

Behavioral Personalization

Reading Distance*
 Near PD*

**Hand-held device required (not included, available upon request)*



m'eyeFit mirror frame clip (included)