



CRIZAL[®] Rx AND STOCK LENS AVAILABILITY CHART

MAY 2018



PLASTIC 1.50				TREXA ⁺			AIRWEAR ⁺			THIN & LITE ⁺								
Clear ¹	Transitions ⁴ **			Clear	Essilor Polarized ^{**}	Transitions ⁴	Clear	Essilor Polarized ^{**}	Transitions ⁴	1.60			1.67 ²			1.74 ³		
	Signature [®] VII	XTRActive [®]	Vantage [®]							Essilor Polarized ^{**}	Signature VII	XTRActive	Vantage	Clear	Essilor Polarized ^{**}	Signature VII	XTRActive	Vantage

PROGRESSIVE LENSES

DEFINITY	DEFINITY [®] 3 PLUS	DEFINITY [®] 3	DEFINITY [®]	DEFINITY SHORT [™]
	ESSILOR PALS	Essilor Ideal [®] Fit [™]	Essilor Ideal [®] Advanced Fit [™]	Essilor Ideal [®]
Essilor Ideal Advanced [™]		Essilor Ideal Advanced Wrap [™]	Essilor Ideal [®] Sport	Essilor Ideal [®] Driving
Essilor Ideal [®] Computer [™]		Accolade [®]	Accolade Freedom [™]	SmallFit [®]
SmallFit [®] Digital		Ovation [®]	Ovation [®] Digital	Essilor Natural [®]
Essilor Natural [®] Digital		Adaptar [®]	Adaptar [®] Digital	Adaptar [®] Short Digital

1 E-SPF[®] of 35 with Crizal Sapphire[®] 360[°] UV and E-SPF of 25 with all other Crizal products made with any lens material other than clear 1.50 plastic.
 2 All 1.67 materials can only be tinted in Gray and Brown.
 3 Thin & Lite 1.74 lenses come systematic with your choice of Crizal[®] Previncia[®], Crizal Sapphire 360[°] UV, Crizal Avancé UV or Crizal Alizé UV.
 4 Transitions products are available in Brown, Gray, and Graphite Green. Please refer to EssilorTransitions.com for details.

■ Available in Essential Blue Series[™] lenses.
 * Crizal Sapphire[®] 360[°] UV available on SFSV Airwear Colors.
 + Crizal Previncia is not available on polarized lenses.
 ▲ Crizal SunShield UV available on polarized and tinted lenses.



PLASTIC 1.50				TREXA ⁺			MID	AIRWEAR ⁺			THIN & LITE ⁺								
Clear ¹	Transitions ⁴ **			Clear	Essilor Polarized ^{**}	Transitions ⁴	1.55 Essilor Polarized ^{**}	Clear	Essilor Polarized ^{**}	Transitions ⁴	1.60			1.67 ²			1.74 ³		
	Signature [®] VII	XTRActive [®]	Vantage [®]								Essilor Polarized ^{**}	Signature VII	XTRActive	Vantage	Clear	Essilor Polarized ^{**}	Signature VII	XTRActive	Vantage

OTHER LENSES

SINGLE VISION	Eyezen [™] + Single Vision	Crizal Kids UV [™]	Crizal [®] Previncia [®] Kids	SV	SV AS	Essilor 360 [°] ™	Essilor eyecode [™] Single Vision	Essilor FIT [™] Single Vision	Essilor Azio [™] Single Vision	Essilor Stylistic [™] WRAP Single Vision	Essilor Sport Solutions [™] 5	
	BIFOCALS	FT28	FT28 AS	FT28 Essilor Bifocal AB [™]	FT35	FT45	ROUND 22					
		TRIFO-CALS	7 x 28	8 x 35								
			MID DISTANCE	Essilor Anti-Fatigue [™]	Essilor Computer [™]	Interview [®]						

1 E-SPF[®] of 35 with Crizal Sapphire[®] 360[°] UV and E-SPF of 25 with all other Crizal products made with any lens material other than clear 1.50 plastic.
 2 All 1.67 materials can only be tinted in Gray and Brown.
 3 Thin & Lite 1.74 lenses come systematic with your choice of Crizal[®] Previncia[®], Crizal Sapphire 360[°] UV, Crizal Avancé UV or Crizal Alizé UV.
 4 Transitions products are available in Brown, Gray, and Graphite Green. Please refer to EssilorTransitions.com for details.
 5 Essilor Sport Solutions has 6 different lens colors available in tint and polarized.

■ Available in Essential Blue Series[™] lenses.
 * Crizal Sapphire[®] 360[°] UV available on SFSV Airwear Colors.
 + Crizal Previncia is not available on polarized lenses.
 ▲ Crizal SunShield UV available on polarized and tinted lenses.



Rx LENS AVAILABILITY
for Eyecare Professionals

Effective MAY 2018



PLASTIC 1.50					TRIVEX ¹				MID	POLY					HI-INDEX							Younger Drivewear [®] 1.50 / Polycarbonate
Clear ¹	Essilor Polarized ⁺⁺	Transitions ⁻³			Clear	Essilor Polarized ⁺⁺	Transitions ⁻³		1.55 Essilor Polarized ⁺⁺	Clear	Essilor Polarized ⁺⁺	Transitions ⁻³			1.60			1.67			1.74 ²	
		Signature [®] VII	XTRActive [®]	Vantage [®]			Essilor Polarized ⁺⁺	Signature VII				XTRActive	Vantage	Clear	Essilor Polarized ⁺⁺	Signature VII	XTRActive	Vantage	Clear ³	Essilor Polarized ⁺⁺		

OTHER LENSES

	Clear ¹	Essilor Polarized ⁺⁺	Signature [®] VII	XTRActive [®]	Vantage [®]	Clear	Essilor Polarized ⁺⁺	Signature VII	XTRActive	Vantage	1.55 Essilor Polarized ⁺⁺	Clear	Essilor Polarized ⁺⁺	Signature VII	XTRActive	Vantage	Clear	Essilor Polarized ⁺⁺	Signature VII	XTRActive	Vantage	Clear ³	Essilor Polarized ⁺⁺	
YOUNGER																								
SV Spherical						•																		•
SV Aspherical						•																		
FT28														•										
Image [®]	•	•	•			•																		•
X-CEL																								
SFSV	•		•	•	•	•		•	•	•														
Asp-SFSV	•					•		•																
D-25	•																							
D-28	•		•	•		•		•																
D-35	•		•			•		•																
7X25	•																							
7X28	•		•	•		•		•																
7X35	•																							
8X35	•		•			•		•																
RD-22	•																							
RD-24			•			•		•																
Blended Rd	•																							
Freedom ID Pal			•			•		•	•	•														

SPECIALTY ITEMS (on UNCOATED LENSES ONLY)

Slaboff, High Adds, Bi-Concave, etc...	•																							
--	---	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--

1 E-SPF[®] of 35 with Crizal Sapphire[®] 360[°] UV and E-SPF of 25 with all other Crizal products made with any lens material other than clear 1.50 plastic.
 2 Thin & Lite 1.74 lenses come systematic with your choice of Crizal[®] Previncia[®], Crizal Sapphire 360[°] UV, Crizal Avancé UV or Crizal Alizé UV.
 3 Transitions products are available in Brown, Gray, and Graphite Green. Please refer to EssilorTransitions.com for details.

• Crizal Sapphire[®] 360[°] UV, not available on polarized or tinted lenses.
 + Crizal Previncia is not available on polarized lenses.
 ▲ Crizal SunShield UV[®] available on polarized and tinted lenses.



Rx LENS AVAILABILITY
for Eyecare Professionals

Effective MAY 2018



FINISHED SINGLE VISION

CRIZAL [®] PREVENCIA [®]			
MATERIAL	INDEX	POWER RANGE	SIZE
Airwear [®] (Spherical)	1.59	+4.00 to -8.00 up to -2.00 cyl -2.00 cyl Available to a -6.00	65/70 mm
Plastic Flat (Spherical)	1.50	+4.00 to -6.00 up to -2.00 cyl -2.00 cyl Available to a -4.00	70 mm

CRIZAL AVANCÉ UV [™]			
MATERIAL	INDEX	POWER RANGE	SIZE
Thin & Lite [®] 1.74 (Aspheric)	1.74	-3.00 to -12.00 up to -2.00 cyl* -2.00 cyl available to a -10.00	70/75/80 mm
Thin & Lite 1.67 (Aspheric)	1.67	+6.00 to -10.00 up to -2.00 cyl 2.00 cyl available to a -10.00	65/70/75/80 mm
Thin & Lite 1.67 Transitions [®] Signature [®] VII Gray (Aspheric)	1.67	-2.00 to -8.00 up to -2.00 cyl -2.00 cyl available to a -8.00	70/75/80 mm
Thin & Lite 1.60 (Aspheric)	1.60	+6.00 to -10.00 up to -2.00 cyl -2.00 cyl available to a -10.00	65/70/75 mm
Airwear [®] (Aspheric) & (Spherical)	1.59	+4.00 to -8.00 up to -2.00 cyl -2.00 cyl available to a -6.00	65/70 mm
Airwear Transitions Signature VII Gray & Graphite Green (Spherical)	1.59	4.00 to -6.00 up to -2.00 cyl -2.00 cyl available to a -4.00	65/70 mm
Plastic Flat (Spherical)	1.50	+4.00 to -6.00 up to -2.00 cyl -2.00 cyl available to a -4.00	70 mm
1.50 Transitions Signature VII Gray (Spherical)	1.50	+4.00 to -4.00 up to -2.00 cyl -2.00 cyl available to a -2.00	65/70 mm

* -10.25 to -12.00 available in sphere only.

FINISHED SINGLE VISION

CRIZAL ALIZÉ UV™			
MATERIAL	INDEX	POWER RANGE	SIZE
Thin & Lite® 1.67 (Spherical)	1.67	-2.00 to -8.00 up to -2.00 cyl -2.00 cyl available to a -6.00	70/75 mm
Thin & Lite 1.67 (Aspheric)	1.67	+6.00 to -12.00 up to -2.00 cyl -2.00 cyl available to a -10.00*	65/70/75/80 mm
Thin & Lite 1.60 (Spherical)	1.60	+6.00 to -10.00 up to -2.00 cyl -2.00 cyl available to a -10.00	65/70/75 mm
Airwear® (Aspheric) & (Spherical)	1.59	+4.00 to -8.00 up to -2.00 cyl -2.00 cyl available to a -6.00	65/70 mm
Airwear® Transitions® Signature® VII Gray (Spherical)	1.59	+4.00 to -6.00 up to -2.00 cyl -2.00 cyl available to a -4.00	65/70 mm
TREXA® (Spherical)	1.53	+4.00 to -6.00 up to -2.00 cyl -2.00 cyl available to a -6.00	65/70 mm
Plastic Flat (Spherical) ¹	1.50	+4.00 to -6.00 up to -2.00 cyl -2.00 cyl available to a -4.00	70 mm
1.50 Transitions® Signature® VII Gray (Spherical)	1.50	+4.00 to -4.00 up to -2.00 cyl -2.00 cyl available to a -2.00	65/70 mm
CRIZAL EASYUV™			
MATERIAL	INDEX	POWER RANGE	SIZE
Thin & Lite 1.67 (Spherical)	1.67	-2.00 to -8.00 up to -2.00 cyl -2.00 cyl available to a -6.00	70/75 mm
Thin & Lite 1.67 (Aspheric)	1.67	-2.00 to -12.00 up to -2.00 cyl -2.00 cyl available to a -10.00*	70/75/80 mm
Airwear (Aspheric) & (Spherical)	1.59	+4.00 to -8.00 up to -2.00 cyl -2.00 cyl available to a -6.00	65/70 mm
Plastic Flat (Spherical)	1.50	+4.00 to -6.00 up to -2.00 cyl -2.00 cyl available to a -4.00	70 mm

* -10.25 to -12.00 available in sphere only.

CRIZAL® NO-GLARE LENSES PROVIDE THE CLEAREST VISION POSSIBLE



Reduce glare

Glare, poor night vision, and reflection are major concerns for patients

- Less reflection during the day and safer nighttime driving



Resist scratches

Scratch resistance is the benefit patients want most in anti-reflective lenses

- Protection against everyday wear and tear, with no scratches from cleaning



Repel water

Patients want lenses that stay clean in outdoor environments

- More convenient to wear outdoors



Resist smudges

Lenses that smudge easily are a major source of frustration

- Fewer smears and fingerprints mean less visual distraction



Repel dust

Patients are bothered by dust, dirt, and lashes clinging to static on the inside of lenses


- *Crizal* repels these particles, so lenses stay cleaner longer



Comprehensive UV protection

UV protection is needed every day, but many patients don't know this

CRIZAL SAPPHIRE® 360° UV LENS 9-STEP PROCESS



E-SPF is an index rating the overall UV protection of a lens. E-SPF was developed by Essilor International and endorsed by third party experts. The E-SPF index relates to lens performance only and excludes direct eye exposure that depends on external factors (wearer's morphology, frame shape, position of wear).

CRIZAL® PRODUCT FAMILY



- More durable than ordinary lenses
- Easier to clean than ordinary lenses

- All the benefits of Crizal Easy UV lenses, PLUS:
- Superior protection against smudges, oil, and water
 - Lenses stay cleaner longer

- All the benefits of Crizal Alizé UV lenses, PLUS:
- 2x more scratch resistance to protect lens investment
 - Repels dust and dirt for clearer vision and less cleaning

- All the benefits of Crizal Avancé UV lenses, PLUS:
- Best-in-class transparency
 - Reduction of reflections from all around the wearer
 - Plus, UV protection has been improved with E-SPF® 35

- Selectively filters at least 20% of Harmful Blue Light*



Applies to the Crizal lens process

CRIZAL[®] NO-GLARE LENSES

Worn and trusted by millions of people worldwide.



Optimal UV protection in a clear lens

The higher the E-SPF[®] Index, the greater the level of UV protection provided by the lens

E-SPF is an index rating the overall UV protection of a lens. *E-SPF* was developed by Essilor International and endorsed by third party experts. The *E-SPF* index relates to lens performance only and excludes direct eye exposure that depends on external factors (wearer's morphology, frame shape, position of wear).

For the most current availability chart, visit [Crizal.com](https://www.crizal.com)



Better Sight.
Better Life.

Transitions™

Crizal[®]

VARILUX[®]

Eyezen™

Xperio^{UV}™