



VARILUX®

LENSES AVAILABILITY for Eyecare Professionals

Effective **NOV 2020**

VARILUX® X DESIGN			LENS ID		VARILUX® X FIT										LENS ID		VARILUX® X 4D			LENS ID	
MATERIAL	PRODUCT	INDEX	G	B	GG	COLORS	G1	B1	G-G	C	R	ESC	POWER RANGE	ADD POWER*	MAX CYLINDER						
PLASTIC	• Clear	1.50											-10.00 to +6.00	0.75 to 4.00D	-6.00						
	Essilor Polarized	1.50	✓	✓									-10.00 to +6.00	0.75 to 4.00D	-6.00						
	Xperio UV™ Polarized	1.50	✓	✓									-10.00 to +6.00	0.75 to 4.00D	-6.00						
	Transitions® Signature® GEN 8™	1.50	✓	✓	✓								-10.00 to +6.00	0.75 to 4.00D	-6.00						
	Transitions® XTRActive®	1.50	✓	✓									-10.00 to +6.00	0.75 to 4.00D	-6.00						
TREXA® (TRIVEX®)	• Clear	1.53											-10.00 to +6.00	0.75 to 4.00D	-6.00						
	Transitions® Signature® GEN 8™	1.53	✓	✓									-10.00 to +6.00	0.75 to 4.00D	-6.00						
AIRWEAR® (POLYCARBONATE)	• Clear	1.59											-10.00 to +6.00	0.75 to 4.00D	-6.00						
	Essilor Polarized	1.59	✓	✓									-10.00 to +6.00	0.75 to 4.00D	-6.00						
	Xperio UV™ Polarized	1.59	✓	✓									-10.00 to +6.00	0.75 to 4.00D	-6.00						
	Transitions® Signature® GEN 8™	1.59	✓	✓	✓	✓							-10.00 to +6.00	0.75 to 4.00D	-6.00						
	Transitions® XTRActive®	1.59	✓SM	✓SM	✓								-10.00 to +6.00	0.75 to 4.00D	-6.00						
THIN & LITE®	• Clear	1.67											-12.00 to +9.00	0.75 to 4.00D	-6.00						
	Transitions® Signature® GEN 8™	1.67	✓	✓	▼	✓							-12.00 to +9.00	0.75 to 4.00D	-6.00						
	Transitions® XTRActive®	1.67	✓SM	✓	✓								-12.00 to +9.00	0.75 to 4.00D	-6.00						
	• Clear**	1.74											-14.00 to +13.00	0.75 to 4.00D	-6.00						
	Transitions® Signature® GEN 8™**	1.74	✓	✓									-14.00 to +13.00	0.75 to 4.00D	-6.00						

MFH=14mm

NOTES: All Varilux digital lenses require a hard coat. Essilor AR stacks are required when No-Glare lens treatment is applied (for example, Crizal® No-Glare lenses).

*ADD Power in steps of 0.25D. At fitting height less than 17mm, ADD power is up to +3.50.

**THIN & LITE 1.74 lenses come systematic with choice of Crizal® Preveincia®, Crizal Sapphire® 360° UV, Crizal Avancé UV™ or Crizal Alizé UV™.

• Available in Essential Blue Series™ lenses.

G: Gray, B: Brown, GG: Graphite Green, Colors: Transitions® Signature® style colors in Sapphire, Amethyst, Amber, or Emerald, G1: Gray 1, B1: Brown 1, G-G: Gray-Green, C: Copper, R: Ruby, ESC: Essilor Sun Colors

✓ **Currently Available**

▼ **Availability date to be announced. Current availability Transitions® Signature® VII Only.**

SM Available in Transitions® XTRActive® Style Mirrors in the following six flash mirror colors: Gold, Silver Shadow, Blue, Green, Pink, and Red.



VARILUX® LENSES AVAILABILITY

for Eyecare Professionals

Effective **NOV 2020**

VARILUX® PHYSIO® W3+														 		VARILUX® PHYSIO® W3+ FIT														 		VARILUX® PHYSIO® W3+ EYECODE™														 	
MATERIAL	PRODUCT	INDEX	G	B	GG	COLORS	G1	B1	G-G	C	R	ESC	POWER RANGE	BASE CURVES	ADD POWER*	CYLINDER																															
PLASTIC	• Clear	1.50											-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50, 4.00D	-6.00																															
	Essilor Polarized	1.50	✓	✓			✓					☼	-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50, 4.00D	-6.00																															
	Xperio UV™ Polarized	1.50	✓	✓			✓					☼	-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50, 4.00D	-6.00																															
	Transitions® Signature® GEN 8™	1.50	✓	✓	✓								-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50, 4.00D	-6.00																															
	Transitions® XTRActive®	1.50	✓										-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50, 4.00D	-6.00																															
TREXA® (TRIVEX®)	• Clear	1.53											-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50, 4.00D	-6.00																															
	Transitions® Signature® GEN 8™	1.53	✓	✓	▼								-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50, 4.00D	-6.00																															
	Transitions® XTRActive®	1.53	✓	✓	✓								-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50, 4.00D	-6.00																															
AIRWEAR® (POLYCARBONATE)	• Clear	1.59											-10.00 to +6.00	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50, 4.00D	-6.00																															
	Essilor Polarized	1.59	✓	✓			✓					☼	-10.00 to +6.00	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50, 4.00D	-6.00																															
	Xperio UV™ Polarized	1.59	✓	✓			✓					☼	-10.00 to +6.00	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50, 4.00D	-6.00																															
	Transitions® Signature® GEN 8™	1.59	✓	✓	▼								-10.00 to +6.00	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50, 4.00D	-6.00																															
	Transitions® XTRActive®	1.59	✓SM	✓SM									-10.00 to +6.00	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50, 4.00D	-6.00																															
THIN & LITE®	Clear	1.60											-12.00 to +8.00	1.75, 3.00, 4.00, 5.50, 7.00, 8.00	0.75 to 3.50, 4.00D	-6.00																															
	Transitions® Signature® GEN 8™	1.60	✓	✓									-12.00 to +8.00	1.75, 3.00, 4.00, 5.50, 7.00, 8.00	0.75 to 3.50, 4.00D	-6.00																															
	Clear	1.67											-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.50, 8.00	0.75 to 3.50, 4.00D	-6.00																															
	Essilor Polarized	1.67	✓	✓									-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.50, 8.00	0.75 to 3.50, 4.00D	-6.00																															
	Xperio UV™ Polarized	1.67	✓	✓									-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.50, 8.00	0.75 to 3.50, 4.00D	-6.00																															
	Transitions® Signature® GEN 8™	1.67	✓	✓	▼								-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.50, 8.00	0.75 to 3.50, 4.00D	-6.00																															
	Transitions® XTRActive®	1.67	✓SM	✓SM	✓								-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.50, 8.00	0.75 to 3.50, 4.00D	-6.00																															
	Clear**	1.74											-18.00 to +13.00	1.75, 2.75, 3.75, 5.25, 6.50, 8.00, 10.00	0.75 to 3.50, 4.00D	-6.00																															
	Transitions® Signature® GEN 8™	1.74	✓										-18.00 to +13.00	0.75, 1.75, 2.75, 3.75, 5.25, 6.50, 8.00, 10.00	0.75 to 3.50, 4.00D	-6.00																															

MFH=14mm

VARILUX® COMFORT MAX														LENS ID		VARILUX® COMFORT MAX FIT														LENS ID	
MATERIAL		PRODUCT		INDEX	G	B	GG	COLORS	G1	B1	G-G	C	R	ESC	POWER RANGE		BASE CURVES		ADD POWER*		CYLINDER										
PLASTIC	• Clear			1.50												-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25		0.75 to 3.50D		-6.00										
	Essilor Polarized			1.50	✓	✓			✓							-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25		0.75 to 3.50D		-6.00										
	Xperio UV™ Polarized			1.50	✓	✓										-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25		0.75 to 3.50D		-6.00										
	Transitions® Signature® GEN 8™			1.50	✓	✓										-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25		0.75 to 3.50D		-6.00										
	Transitions® XTRActive®			1.50	✓											-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25		0.75 to 3.50D		-6.00										
TREXA® (TRIVEX®)	• Clear			1.53												-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25		0.75 to 3.50D		-6.00										
	Transitions® Signature® GEN 8™			1.53	✓	✓										-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25		0.75 to 3.50D		-6.00										
	Transitions® XTRActive®			1.53	✓											-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25		0.75 to 3.50D		-6.00										
AIRWEAR® (POLYCARBONATE)	• Clear			1.59												-10.00 to +6.00	1.75, 3.00, 4.00, 5.50, 7.25		0.75 to 3.50D		-6.00										
	Essilor Polarized			1.59	✓	✓			✓							-10.00 to +6.00	1.75, 3.00, 4.00, 5.50, 7.25		0.75 to 3.50D		-6.00										
	Xperio UV™ Polarized			1.59	✓	✓										-10.00 to +6.00	1.75, 3.00, 4.00, 5.50, 7.25		0.75 to 3.50D		-6.00										
	Transitions® Signature® GEN 8™			1.59	✓	✓			✓							-10.00 to +6.00	1.75, 3.00, 4.00, 5.50, 7.25		0.75 to 3.50D		-6.00										
	Transitions® XTRActive®			1.59	✓SM	✓SM										-10.00 to +6.00	1.75, 3.00, 4.00, 5.50, 7.25		0.75 to 3.50D		-6.00										
THIN & LITE®	Clear			1.60				✓								-12.00 to +8.00	1.75, 3.00, 4.00, 5.50, 7.00, 8.00		0.75 to 3.50D		-6.00										
	Transitions® Signature® GEN 8™			1.60	✓	✓										-12.00 to +8.00	1.75, 3.00, 4.00, 5.50, 7.00, 8.00		0.75 to 3.50D		-6.00										
	• Clear			1.67												-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.50, 8.00		0.75 to 3.50D		-6.00										
	Essilor Polarized			1.67	✓	✓										-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.50, 8.00		0.75 to 3.50D		-6.00										
	Xperio UV™ Polarized			1.67	✓	✓										-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.50, 8.00		0.75 to 3.50D		-6.00										
	Transitions® Signature® GEN 8™			1.67	✓	✓			✓							-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.50, 8.00		0.75 to 3.50D		-6.00										
	Transitions® XTRActive®			1.67	✓SM											-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.50, 8.00		0.75 to 3.00D		-6.00										

MFH=14mm

NOTES: All Varilux digital lenses require a hard coat. Essilor AR stacks are required when No-Glare lens treatment is applied (for example, Crizal® No-Glare lenses).

*ADD Power in steps of 0.25D. At fitting height less than 17mm, ADD power is up to +3.00.

**THIN & LITE 1.74 lenses come systematic with choice of Crizal® Prevencia®, Crizal Sapphire® 360° UV, Crizal Avancé UV™ or Crizal Alizé UV™.

• Available in Essential Blue Series™ lenses.

G: Gray, B: Brown, GG: Graphite Green, Colors: Transitions® Signature® style colors in Sapphire, Amethyst, Amber, or Emerald, G1: Gray 1, B1: Brown 1, G-G: Gray-Green, C: Copper, R: Ruby, ESC: Essilor Sun Colors

✓ Currently Available

▼ Availability date to be announced. Current availability Transitions® Signature® VII Only.

SM Available in Transitions® XTRActive® Style Mirrors in the following six flash mirror colors: Gold, Silver Shadow, Blue, Green, Pink, and Red.

☼ Essilor Sun Colors (solid or gradient: gray mist, sienna brown, forest gray, grape, plum, emerald blue)



VARILUX® LENSES AVAILABILITY

for Eyecare Professionals

Effective **NOV 2020**

VARILUX PHYSIO DRx™													LENS ID H → 0			
MATERIAL	PRODUCT	INDEX	G	B	GG	COLORS	G1	B1	G-G	C	R	ESC	POWER RANGE	BASE CURVES	ADD POWER*	CYLINDER
PLASTIC	• Clear	1.50											-10.00 to +7.50	2.00, 3.50, 4.25, 6.25, 7.25, 8.75	0.75 to 3.50, 4.00D	-6.00
	Essilor Polarized	1.50	✓	✓			✓	✓	✓			☼	-10.00 to +7.50	2.00, 3.25, 4.25, 6.25, 7.25, 8.25	0.75 to 3.50, 4.00D	-6.00
	Xperio UV™ Polarized	1.50	✓	✓			✓	✓	✓			☼	-10.00 to +7.50	2.00, 3.25, 4.25, 6.25, 7.25, 8.25	0.75 to 3.50, 4.00D	-6.00
	Transitions® Signature® GEN 8™	1.50	✓	✓	✓	✓							-10.00 to +7.50	2.00, 3.50, 4.25, 6.25, 7.25, 8.75	0.75 to 3.50, 4.00D	-6.00
	Transitions® XTRActive®	1.50	✓	✓	✓								-10.00 to +6.00	2.25, 4.25, 5.25, 8.25	0.75 to 3.50, 4.00D	-6.00
	Transitions® Vantage®	1.50	✓											-10.00 to +6.00	2.25, 4.25, 5.25, 8.25	0.75 to 3.50, 4.00D
TREXA® (TRIVEX®)	• Clear	1.53											-10.00 to +7.50	2.00, 3.50, 4.25, 6.25, 7.25, 8.75	0.75 to 3.50, 4.00D	-6.00
	Transitions® Signature® GEN 8™	1.53	✓	✓	▼								-10.00 to +7.50	2.00, 3.50, 4.25, 6.25, 7.25, 8.75	0.75 to 3.50, 4.00D	-6.00
	Transitions® XTRActive®	1.53	✓	✓	✓								-10.00 to +7.50	2.00, 3.50, 4.25, 6.25, 7.25, 8.75	0.75 to 3.50, 4.00D	-6.00
AIRWEAR® (POLYCARBONATE)	• Clear	1.59											-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	0.75 to 3.50, 4.00D	-6.00
	Essilor Polarized	1.59	✓	✓			✓	✓	✓	✓	✓	☼	-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	0.75 to 3.50, 4.00D	-6.00
	Xperio UV™ Polarized	1.59	✓	✓			✓	✓	✓	✓	✓	☼	-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	0.75 to 3.50, 4.00D	-6.00
	Transitions® Signature® GEN 8™	1.59	✓	✓	✓	✓							-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	0.75 to 3.50, 4.00D	-6.00
	Transitions® XTRActive®	1.59	✓SM	✓SM	✓								-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	0.75 to 3.50, 4.00D	-6.00
	Transitions® Vantage®	1.59	✓										-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	0.75 to 3.50, 4.00D	-6.00
THIN & LITE®	Clear	1.60											-12.00 to +8.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.50	0.75 to 3.50, 4.00D	-6.00
	Transitions® Signature® GEN 8™	1.60	✓	✓									-12.00 to +8.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.50	0.75 to 3.50, 4.00D	-6.00
	Clear	1.67											-14.00 to +9.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.00	0.75 to 3.50, 4.00D	-6.00
	Essilor Polarized	1.67	✓	✓					✓				-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.75, 8.00	0.75 to 3.50, 4.00D	-6.00
	Xperio UV™ Polarized	1.67	✓	✓					✓				-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.75, 8.00	0.75 to 3.50, 4.00D	-6.00
	Transitions® Signature® GEN 8™	1.67	✓	✓	✓	✓							-14.00 to +9.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.00	0.75 to 3.50, 4.00D	-6.00
	Transitions® XTRActive®	1.67	✓SM	✓SM	✓								-14.00 to +9.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.00	0.75 to 3.50, 4.00D	-6.00
	Clear**	1.74											-20.00 to +13.00	0.75, 1.75, 2.75, 3.75, 5.25, 6.50, 8.00, 10.00	0.75 to 3.50, 4.00D	-6.00
	Transitions® Signature® GEN 8™	1.74	✓										-20.00 to +13.00	0.75, 1.75, 2.75, 3.75, 5.25, 6.50, 8.00, 10.00	0.75 to 3.50, 4.00D	-6.00

MFH=17mm

NOTES: All Varilux digital lenses require a hard coat. Essilor AR stacks are required when No-Glare lens treatment is applied (for example, Crizal® No-Glare lenses).

**THIN & LITE 1.74 lenses come systematic with choice of Crizal® Prevenir®, Crizal Sapphire® 360° UV, Crizal Avancé UV™ or Crizal Alizé UV™.

***For the complete list of Transitions color availability please refer to the Transitions Availability Chart.

For complete Sunwear availability please refer to the Essilor Sun Portfolio Availability Chart.

• Available in Essential Blue Series™ lenses.

G: Gray, B: Brown, GG: Graphite Green, Colors: Transitions® Signature® style colors in Sapphire, Amethyst, Amber, or Emerald, G1: Gray 1, B1: Brown 1, G-G: Gray-Green, C: Copper, R: Ruby, ESC: Essilor Sun Colors

✓ **Currently Available**

▼ **Availability date to be announced. Current availability Transitions® Signature® VII Only.**

SM **Available in Transitions® XTRActive® Style Mirrors in the following six flash mirror colors: Gold, Silver Shadow, Blue, Green, Pink, and Red.**

☼ **Essilor Sun Colors (solid or gradient: gray mist, sienna brown, forest gray, grape, plum, emerald blue)**

VARILUX PHYSIO SHORT DRx™

Rx availability for Varilux Physio Short DRx matches Varilux Physio DRx, except ADD power range is 0.75 to 3.00D across all materials.

MFH=14mm



VARILUX® LENSES AVAILABILITY

for Eyecare Professionals

Effective **NOV 2020**

VARILUX COMFORT DRx™							
MATERIAL	PRODUCT	INDEX	COLOR	POWER RANGE	BASE CURVES	ADD POWER	CYLINDER
PLASTIC	● 1.50 Index Clear	1.50		-10.00 to +6.00	2.00, 3.50, 4.25, 6.25, 7.25	0.75 to 3.50D	-6.00
	1.50 Index Essilor Polarized	1.50	Gray, Brown, Gray-Green	-10.00 to +6.00	2.00, 3.25, 4.25, 6.25, 8.25	0.75 to 3.50D	-4.00
	1.50 Index Xperio UV™ Polarized	1.50	Gray, Brown, Gray-Green	-10.00 to +6.00	2.00, 3.25, 4.25, 6.25, 8.25	0.75 to 3.50D	-4.00
	1.50 Index Transitions® Signature® VII	1.50	Graphite Green, Gray, Brown	-10.00 to +6.00	2.00, 3.50, 4.25, 6.25, 7.25	0.75 to 3.50D	-6.00
	1.50 Index Transitions® XTRActive®	1.50	Graphite Green, Gray, Brown	-10.00 to +6.00	2.25, 4.25, 5.25, 8.25	0.75 to 3.50D	-6.00
	1.50 Index Transitions® Vantage®	1.50	Gray	-10.00 to +6.00	2.25, 4.25, 5.25, 8.25	0.75 to 3.50D	-6.00
TREXA® (TRIVEX®)	● Clear	1.53		-10.00 to +6.00	2.00, 3.50, 4.25, 6.25, 7.25	0.75 to 3.50D	-6.00
	Transitions® Signature® VII	1.53	Graphite Green, Gray, Brown	-10.00 to +6.00	2.00, 3.50, 4.25, 6.25, 7.25	0.75 to 3.50D	-6.00
	Transitions® XTRActive®	1.53	Graphite Green, Gray, Brown	-10.00 to +6.00	2.00, 3.50, 4.25, 6.25, 7.25	0.75 to 3.50D	-6.00
AIRWEAR® (POLYCARBONATE)	● Clear	1.59		-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.50, 8.50	0.75 to 3.50D	-6.00
	Essilor Polarized	1.59	Gray, Brown, Gray-Green	-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.50, 8.50	0.75 to 3.50D	-6.00
	Xperio UV™ Polarized	1.59	Gray, Brown, Gray-Green, Copper, Ruby	-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.50, 8.50	0.75 to 3.50D	-6.00
	Transitions® Signature® VII	1.59	Graphite Green, Gray, Brown	-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.50, 8.50	0.75 to 3.50D	-6.00
	Transitions® XTRActive®	1.59	Graphite Green, Gray, Brown	-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.50, 8.50	0.75 to 3.50D	-6.00
	Transitions® Vantage®	1.59	Gray	-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.50, 8.50	0.75 to 3.50D	-6.00
THIN&LITE®	● 1.60 Clear	1.60		-12.00 to +8.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.50	0.75 to 3.50D	-6.00
	1.60 Transitions® Signature® VII	1.60	Gray, Brown	-12.00 to +8.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.50	0.75 to 3.50D	-6.00
	● 1.67 Clear	1.67		-14.00 to +9.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.00	0.75 to 3.50D	-6.00
	1.67 Essilor Polarized	1.67	Gray, Brown, Gray-Green	-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.75, 8.00	0.75 to 3.50D	-6.00
	1.67 Xperio UV™ Polarized	1.67	Gray, Brown, Gray-Green	-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.75, 8.00	0.75 to 3.50D	-6.00
	1.67 Transitions® Signature® VII	1.67	Graphite Green, Gray, Brown	-14.00 to +9.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.00	0.75 to 3.50D	-6.00
	1.67 Transitions® XTRActive®	1.67	Graphite Green, Gray, Brown	-14.00 to +9.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.00	0.75 to 3.50D	-6.00
	● 1.74 Clear*	1.74		-20.00 to +13.00	0.75, 1.75, 2.75, 3.75, 5.25, 6.50, 8.00, 10.00	0.75 to 3.50D	-6.00
	1.74 Transitions® Signature® VII	1.74	Gray	-20.00 to +13.00	0.75, 1.75, 2.75, 3.75, 5.25, 6.50, 8.00, 10.00	0.75 to 3.50D	-6.00

MFH=17mm

NOTES: All Varilux digital lenses require a hard coat. Essilor AR stacks are required when No-Glare coating is applied (for example, Crizal® No-Glare lenses).

*Thin&Lite 1.74 lenses come systematic with choice of Crizal® Previncia®, Crizal Sapphire® 360° UV, Crizal Avancé UV™ or Crizal Alizé UV™.

Please see Xperio UV™ availability chart for full range of Xperio UV colors including Essilor® Sunwear Colors™.

Trivex is a registered trademark of PPG Industries Ohio, Inc.

- Available in Essential Blue Series™ lenses.

VARILUX COMFORT SHORT DRx™		LENS ID	

Rx availability for Varilux Comfort Short DRx matches Varilux Comfort DRx, except ADD power range is 0.75 to 3.00D across all materials.

MFH=14mm

Transitions™

Crizal®

VARILUX®

Eyezen™

Xperio UV™

